

Pertex®



| REF | |
|----------|----------|
| 00801 | 00811-EX |
| 00801-EX | 00840-05 |
| 00811 | 00840-EX |

INTENDED PURPOSE

Pertex® is a mounting medium intended for use in in vitro diagnostics for the mounting of histological and cytological human tissue samples on slides. Mounting with Pertex® enables pathological diagnosis and is intended for use by trained personnel in pathology laboratories.

USE

Pertex® is a xylene-based and optically clear mounting agent used in in vitro diagnostics for mounting histological and cytological human tissue samples on slides. Mounting with Pertex® enables examination of the tissue samples under an optical microscope for diagnosis and ensures, through the airtight cover, that the tissue sample retains its colour, shape and structure during long-term storage.

To visualise the different structures and cellular details of tissue or cell samples, thin sections or smears of these are stained on a slide with histological or cytological staining solutions. The staining solutions are usually water-based, and the tissue is thus filled with water after staining. In the final stages of the staining process, the slide containing the tissue is passed through baths of increasing concentration of alcohol to progressively dehydrate the tissue. The final step is a bath of organic solvents.

Pertex® consists of synthetic resin dissolved in xylene. Pertex® is dripped onto the slide after the clearing step, the cover slip is mounted and the xylene evaporates, hardening and forming a hard, airtight, clear film that bonds the cover slip to the slide and protects the tissue. The refractive index of Pertex® (1.492) is close to that of glass, and the adhesive will therefore not be visible or interfere with examination of the preparations under the microscope.

Pertex® is fast drying, enabling immediate examination of the preparations. It also has low viscosity, which allows it to flow easily over the glass and prevents bubbles.

SPECIFICATION

Composition

| | |
|--------------|-------------|
| Xylene | 30-65 % |
| Ethylbenzene | 0-20 % |
| Toluene | <0,5 % |
| Resin | Proprietary |

Properties

| | |
|------------------|-------------|
| Refractive index | 1,492 |
| Viscosity | 650 ± 50 cP |

INSTRUCTIONS FOR USE

Storage and shelf life

Store in tightly sealed original packaging in a dry and cool place. Store in a well-ventilated place. Protect from heat, sparks and naked flames. Protect from sunlight. Keep away from sources of ignition – smoking prohibited. Store separately from oxidising agents.

The product has a shelf life of 3 years from the date of manufacture in unopened packaging, tubes -shelf life of 5 years. When the packaging is opened, the xylene evaporates, making the solution more viscous over time. The expiry date is printed on the package label.

Warnings/precautions for safe handling

Classification and labelling information in accordance with Regulation (EC) No 1272/2008 (CLP) can be found on the product label and/or safety data sheet.

Wear tight-fitting safety goggles or a face shield. Protective gloves must be worn. If there is a risk of skin contact, appropriate protective clothing should be worn. Use respiratory protection in case of insufficient ventilation.

Ensure good ventilation. Avoid inhaling fumes. Work with point extraction or in a fume cupboard. Avoid contact with skin and eyes. Prevent sparking due to static electricity. Protect from heat, sparks and naked flames. Eliminate all sources of ignition. Fumes can collect at floors and in low-lying areas. Use non-sparking tools and explosion-proof equipment.

Waste management

The product is classified as hazardous waste. Present for destruction according to local regulations. See safety data sheet for more information.

Sample material

The various types of samples that can be analysed are histological and cytological tissue samples.

Preparations

Pertex® is ready to use and does not need to be diluted or modified.

Instructions

1. After clearing in xylene, it is possible to mount immediately.
2. If isopropanol or HistoLab Clear has been used as a clearing agent, allow the glasses to evaporate slightly before mounting.
3. Place the slide with the section horizontal.
4. Drip Pertex® onto the slide. Pertex® flows easily over the glass.
5. Place the cover slip on the slide from one side and carefully fold it over the section to the other side to minimise air bubbles.
6. Leave to dry in a fume cupboard.
7. Allow to dry for about a day before archiving.

The cover slip and the mounting medium can be removed by placing the glass in xylene and leaving it preferably overnight.



ADDITIONAL INFORMATION

Use equipment and reagents suitable for in vitro diagnostics

All serious incidents that have occurred in connection with the product must be reported to the manufacturer and the competent authority.

VERSION HISTORY

Changes since previous version:

Inserted information about intended purpose.

Information about storage.

Information about precautions for safe handling.

Addendum on what to do in case of serious incidents.

SOURCES

Histotechnology A Self Instructional Text, Freida L. Carson

Version / datum: v2-2022-05

| | | | | | |
|---|---|---|---|---|---|
|  |  |  |  |  |  |
| Catalog No. | Batch No. | In Vitro Diagnostic Use | Use By | Consult Instructions for Use | Manufacturer |

Histolab Products AB

Södra Långebergsgatan 36
SE-436 32 Askim, Sweden
P: +46 31 709 30 30
mail@histolab.se
www.histolab.se

